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
鎳鈦旋轉器械的修型概念（Ni-Ti rotary instrument shaping concept）

- ◇ Single length system
- ◇ Crown down system

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
	Color	Tip Size	Taper (%)	Taper Indicator	Cross-Section	Length (mm)	Suggested Angle
E-FLEX REC	●	20	4			21/25/31	FWD: 30° REV: 150°
	●	25	4				
	●	35	4				
	●	40	4				
	●	20	6				
	●	25	6				
	●	35	6				
	●	40	6				

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直的根管 (Straight)	step back / step down / crown down	Single length
彎曲根管 (curved)	step down / crown down Anticurvature filling	Crown down/ Anticurvature filling
牙髓炎 (Normal/Pulpitis)	step back / step down / crown down	Single length
牙髓壞死 (necrosis)	step down / crown down	Crown down
橢圓型根管 (oval)	Circumferential filling envelope of motion Shaping as 2 canals	Anatomic shaping Shaping as 2 canals
C型根管系統 (C-shape)	Circumferential filling envelope of motion	Sufficient cleaning system

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1. Schilder, H., *Cleaning and shaping the root canal*. Dental clinics of north America, 1974. 18 (2) : p. 269-296.
2. Abou-Rass, M., A.L. Frank, and D.H. Glick, *The anticurvature filing method to prepare the curved root canal*. Journal of the American Dental Association (1939), 1980. 101 (5) : p. 792-794.
3. Roane, J.B., C.L. Sabala, and M.G. Duncanson Jr, *The "balanced force" concept for instrumentation of curved canals*. Journal of endodontics, 1985. 11 (5) : p. 203-211.
4. Lim, S. and C. Stock, *The risk of perforation in the curved canal: anticurvature filing compared with the stepback technique*. International endodontic journal, 1987. 20 (1) : p. 33-39.
5. Donald, C.Y., A. Tam, and F.H. Schilder, *Patency and envelope of motion—Two essential procedures for cleaning and shaping the root canal systems*. General dentistry, 2009. 57 (6) : p. 616-621.
6. Kerekes, K. and L. Tronstad, *Long-term results of endodontic treatment performed with a standardized technique*. Journal of endodontics, 1979. 5 (3) : p. 83-90.
7. Walton, R.E., *Histologic evaluation of different methods of enlarging the pulp canal space*. Journal of endodontics, 1976. 2 (10) : p. 304-311.
8. Marshall, F., *A crown-down pressureless preparation root canal enlargement technique, technique manual*. Portland, Ore, 1980.
9. Morgan, L.F. and S. Montgomery, *An evaluation of the crown-down pressureless technique*. Journal of Endodontics, 1984. 10 (10) : p. 491-498.
10. Walia, H., W.A. Brantley, and H. Gerstein, *An initial investigation of the bending and*

- torsional properties of Nitinol root canal files*. Journal of endodontics, 1988. 14 (7) : p. 346-351.
11. Laurindo, F.V. and J.A.P. de Figueiredo, *Reciprocating versus Rotary instruments: a review*. Revista Odonto Ciencia, 2016. 31 (3) : p. 135-139.
12. Bürklein, S. and E. Schäfer, *Apically extruded debris with reciprocating single-file and full-sequence rotary instrumentation systems*. Journal of endodontics, 2012. 38 (6) : p. 850-852.
13. Pedullà, E., et al., *Influence of continuous or reciprocating motion on cyclic fatigue resistance of 4 different nickel-titanium rotary instruments*. Journal of endodontics, 2013. 39 (2) : p. 258-261.
14. Hamze, F., K. Honardar, and K. Nazarimoghadam, *Comparison of two canal preparation techniques using mtwo rotary instruments*. Iranian endodontic journal, 2011. 6 (4) : p. 150.

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E-FLEX S

◆ 3 Shaping files for a deep shape.

◆ 3 Finishing files for navigating challenging curves in the apical region.



Safety Non-cutting Tip

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• Usage Protocol

綜合版 : SX、S1、S2、F1、F2、F3

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Taper/
Cross
Section

Length

Speed/
Torque

#16 #18 #20 #25 #30

Variable Taper

19mm 21/25/31mm

350rpm 2 N-cm

E-FLEX BLUE


Blue Heat Treatment provides superior strength while maintaining flexibility to negotiate all canals

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• Usage Protocol

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4支裝	#17	#15	#20	#25	#30	#35	#40
Taper/ Cross Section	8%	4%	6%				
Length	19mm	21/25/31mm					
Speed/ Torque	300rpm 3N-cm	350rpm 2.5N-cm					

E-FLEX REC


One file reciprocation for incredible efficiency

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• Usage Protocol

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4支裝	#20	#25	#35	#40
Taper/ Cross Section	4%	6%		
Length	21/25/31mm			
Speed/ Torque	FWD:30° REV:150°			

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Young Dental Co., Ltd.